

DNMT3B

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 96kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human DNMT3B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ICF; ICF1; M.HsaIIIB;

Catalog #: A11079

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1789

Isotype: IgG

Swiss Prot: Q9UBC3

Purity: Affinity purification

For research use only.

Background:

CpG methylation is an epigenetic modification that is important for embryonic development, imprinting, and X-chromosome inactivation. Studies in mice have demonstrated that DNA methylation is required for mammalian development. This gene encodes a DNA methyltransferase which is thought to function in de novo methylation, rather than maintenance methylation. The protein localizes primarily to the nucleus and its expression is developmentally regulated. Mutations in this gene cause the immunodeficiency-centromeric instability-facial anomalies (ICF) syndrome. Eight alternatively spliced transcript variants have been described. The full length sequences of variants 4 and 5 have not been determined. [provided by RefSeq, May 2011]

To place an order, please [Click HERE](#).