

PIK3R2

Reactivity: Human Rat

Tested applications: WB IHC IF IP FC

Observed MW: Refer to figures

Immunogen:

Recombinant protein

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A11086

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 5296

Isotype: IgG

Swiss Prot: O00459

Purity: Affinity purification

For research use only.