

CYR61

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:1000 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 42kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CYR61

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

a

Synonym:

CYR61;IGFBP10;IGF binding protein 10; IBP 10; CCN family member 1; CCN1

Catalog #: A1111

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3491

Isotype: IgG

Swiss Prot: O00622

Purity: Affinity purification

For research use only.

Background:

Cyr61 is a secreted heparin binding protein, encoded by a growth factor-inducible immediate-early gene, that associates with the extracellular matrix and connective tissue. Cyr61 is a member of a distinct family of angiogenic and vasculogenic regulators designated CCN proteins, which includes connective tissue growth factor (CTGF) and the mouse Cyr61 homolog, Fisp12. As an angiogenic inducer, Cyr61 binds to the cell surface receptor Integrin α V/ β 3, where it then stimulates cell adhesion and migration and promotes DNA synthesis of human vascular endothelial cells. Expression of Cyr61 is elevated during vessel growth, wound healing and chondrocyte differentiation. Cyr61 is also detected in a wide variety of tumors as it induces tumor growth and functions as a marker of tumor progression.

To place an order, please [Click HERE](#).