

ALAS1

Reactivity: Human

Tested applications: WB ICC IF

Recommended Dilution: WB 1:500 - 1:2000 ICC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 71 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ALAS1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ALAS; MIG4; ALAS3; ALASH; ALAS-H;

Catalog #: A11120

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 211

Isotype: IgG

Swiss Prot: P13196

Purity: Affinity purification

For research use only.

Background:

This gene encodes the mitochondrial enzyme which catalyzes the rate-limiting step in heme (iron-protoporphyrin) biosynthesis. The enzyme encoded by this gene is the housekeeping enzyme; a separate gene encodes a form of the enzyme that is specific for erythroid tissue. The level of the mature encoded protein is regulated by heme: high levels of heme down-regulate the mature enzyme in mitochondria while low heme levels up-regulate. A pseudogene of this gene is located on chromosome 12. Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).