

CRADD

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 23kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CRADD

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

CRADD;MGC9163;RAIDD ;

Catalog #: A1124

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8738

Isotype: IgG

Swiss Prot: P78560

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a death domain (CARD/DD)-containing protein and has been shown to induce cell apoptosis. Through its CARD domain, this protein interacts with, and thus recruits, caspase 2/ICH1 to the cell death signal transduction complex that includes tumor necrosis factor receptor 1 (TNFR1A), RIPK1/RIP kinase, and numbers of other CARD domain-containing proteins. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).