

VTN

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 54/75 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human VTN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VN; V75; VNT;

Catalog #: A11245

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 7448

Isotype: IgG

Swiss Prot: P04004

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond.

To place an order, please [Click HERE](#).