

RAB9A

Reactivity:Human Mouse Rat

Tested applications:WB ICC IF

Recommended Dilution:WB 1:500 - 1:2000 ICC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:23 kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RAB9A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAB9;

To place an order, please [Click HERE](#).

Catalog #:A11290

Antibody Type:

Monoclonal Antibody

Species:Rabbit

Gene ID:9367

Isotype:IgG

Swiss Prot:P51151

Purity:Affinity purification

For research use only.