

KIF23

Reactivity: Human Mouse Rat

Tested applications: WB ICC IF

Recommended Dilution: WB 1:500 - 1:2000 ICC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 98 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KIF23

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CHO1; KNSL5; MKLP1; MKLP-1;

Catalog #: A11293

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 9493

Isotype: IgG

Swiss Prot: Q02241

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of kinesin-like protein family. This family includes microtubule-dependent molecular motors that transport organelles within cells and move chromosomes during cell division. This protein has been shown to cross-bridge antiparallel microtubules and drive microtubule movement in vitro. Alternate splicing of this gene results in multiple transcript variants.

To place an order, please [Click HERE](#).