

TOMM20

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF IP FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100
IP 1:20 - 1:50 FC 1:20 - 1:50

Calculated MW: 16 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TOMM20

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAS20; MOM19; TOM20;

To place an order, please [Click HERE](#).

Catalog #: A11308

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 9804

Isotype: IgG

Swiss Prot: Q15388

Purity: Affinity purification

For research use only.