

ARPC5

Reactivity: Human Mouse Rat

Tested applications: WB ICC IF IP FC

Recommended Dilution: WB 1:500 - 1:2000 ICC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:20 - 1:50
FC 1:20 - 1:50

Calculated MW: 16 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ARPC5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARC16; p16-Arc; dJ127C7.3;

Catalog #: A11317

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 10092

Isotype: IgG

Swiss Prot: O15511

Purity: Affinity purification

For research use only.

Background:

This gene encodes one of seven subunits of the human Arp2/3 protein complex. The Arp2/3 protein complex has been implicated in the control of actin polymerization in cells and has been conserved through evolution. The exact role of the protein encoded by this gene, the p16 subunit, has yet to be determined. Alternatively spliced transcript variants encoding different isoforms have been observed for this gene.

To place an order, please [Click HERE](#).