

SIRT1

Reactivity: Human Mouse

Tested applications: WB IHC IF

Observed MW: Refer to figures

Immunogen:

Synthetic Peptide

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A11362

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 23411

Isotype: IgG

Swiss Prot: Q96EB6

Purity: Affinity purification

For research use only.