

RBFOX3

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:200
FC 1:20 - 1:50

Calculated MW: 34 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RBFOX3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FOX3; NEUN; FOX-3; HRNBP3;

To place an order, please [Click HERE](#).

Catalog #: A11449

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 146713

Isotype: IgG

Swiss Prot: A6NFN3

Purity: Affinity purification

For research use only.