

MAD1L1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 83kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MAD1L1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAD1L1;HsMAD1;MAD1;PIG9;TP53I9;TXBP181 ;

Catalog #: A1153

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8379

Isotype: IgG

Swiss Prot: Q9Y6D9

Purity: Affinity purification

For research use only.

Background:

MAD1L1 is a component of the mitotic spindle-assembly checkpoint that prevents the onset of anaphase until all chromosome are properly aligned at the metaphase plate. MAD1L1 functions as a homodimer and interacts with MAD2L1. MAD1L1 may play a role in cell cycle control and tumor suppression. Three transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).