

OSTM1

Reactivity: Human Rat

Tested applications: WB

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human OSTM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A11698

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 28962

Isotype: IgG

Swiss Prot: Q86WC4

Purity: Affinity purification

For research use only.