

CYP1B1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 61kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CYP1B1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CP1B; GLC3A; P4501B1; CYP1B1;

Catalog #: A1377

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1545

Isotype: IgG

Swiss Prot: Q16678

Purity: Affinity purification

For research use only.

Background:

Cytochrome P450 1B1 (CYP1B1) is a key enzyme involved in the production of potentially carcinogenic estrogen metabolites and the activation of environmental carcinogens and is the predominant member of the CYP1 family expressed in normal breast tissue and breast cancer. Estrogen has been proposed to trigger breast cancer development via an initiating mechanism involving its metabolite, catechol estrogen (CE). CYP1B1 catalyzes the conversion of 17-beta-estradiol to the catechol estrogen metabolites 2-OH-E2 and 4-OH-E2 which have both been postulated to be involved in mammary carcinogenesis. Genetic polymorphisms in CYP1B1 may play an important role in human prostate carcinogenesis as well. Polymorphism of the CYP1B1 gene at codon 432 (Val-->Leu) is associated with a change in catalytic function.

To place an order, please [Click HERE](#).