

IL12B

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 41kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human IL12B

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CLMF; NKSF; CLMF2; NKSF2; IL-12B;

Catalog #: A1394

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3593

Isotype: IgG

Swiss Prot: P29460

Purity: Affinity purification

For research use only.

Background:

IL-12 beta is a subunit of interleukin 12, a cytokine that acts on T and natural killer cells, and has a broad array of biological activities. Overexpression of this protein was observed in the central nervous system of patients with multiple sclerosis (MS), suggesting a role of this cytokine in the pathogenesis of the disease. The promoter polymorphism of IL-12 has been reported to be associated with the severity of atopic and non-atopic asthma in children (1).

To place an order, please [Click HERE](#).