

PRDX4

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 31kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PRDX4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AOE37-2; PRX-4;

Catalog #: A1486

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10549

Isotype: IgG

Swiss Prot: Q13162

Purity: Affinity purification

For research use only.

Background:

PRDX4(Peroxiredoxin-4) is also named as AOE37-2, Prx-IV and belongs to the AhpC/TSA family. PRDX4 is associated with acrosome formation during rat spermatogenesis and has a protective role in the male reproductive tract, because the phenotype of mice lacking this isoform includes testicular atrophy and increased sperm DNA damage(PMID:20864641). It is initially synthesized as a membrane-binding 31-kDa protein and processed into a 27-kDa secretory form and is discarded with the residual bodies(PMID:19208552). PRDX4 is a pentamer of dimers(PMID:23025503). This antibody is specific to PRDX4.

To place an order, please [Click HERE](#).