

CDK9

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:1000 - 1:3000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 43kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CDK9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TAK; C-2k; CTK1; CDC2L4; PITALRE

Catalog #: A1564

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1025

Isotype: IgG

Swiss Prot: P50750

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the cyclin-dependent protein kinase (CDK) family. CDK family members are highly similar to the gene products of *S. cerevisiae* cdc28, and *S. pombe* cdc2, and known as important cell cycle regulators. This kinase was found to be a component of the multiprotein complex TAK/P-TEFb, which is an elongation factor for RNA polymerase II-directed transcription and functions by phosphorylating the C-terminal domain of the largest subunit of RNA polymerase II. This protein forms a complex with and is regulated by its regulatory subunit cyclin T or cyclin K. HIV-1 Tat protein was found to interact with this protein and cyclin T, which suggested a possible involvement of this protein in AIDS.

To place an order, please [Click HERE](#).