

UGT2B7

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 60kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human UGT2B7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UGT2B9; UDPGTH2; UDPGT2B7; UDPGT 2B9;

Catalog #: A1693

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7364

Isotype: IgG

Swiss Prot: P16662

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the UDP-glycosyltransferase (UGT) family. UGTs serve a major role in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds. This protein is localized in the microsome membrane, and has unique specificity for 3,4-catechol estrogens and estriol, suggesting that it may play an important role in regulating the level and activity of these potent estrogen metabolites.

To place an order, please [Click HERE](#).