

MAD2L1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 24kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MAD2L1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAD2; HSMAD2

Catalog #: A1699

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4085

Isotype: IgG

Swiss Prot: Q13257

Purity: Affinity purification

For research use only.

Background:

MAD2L1 is a component of the mitotic spindle assembly checkpoint that prevents the onset of anaphase until all chromosomes are properly aligned at the metaphase plate. MAD2L1 is related to the MAD2L2 gene located on chromosome 1. A MAD2 pseudogene has been mapped to chromosome 14.

To place an order, please [Click HERE](#).