

IGFBP5

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 31kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human IGFBP5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IGFBP5; IBP5;

Catalog #: A1720

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3488

Isotype: IgG

Swiss Prot: P24593

Purity: Affinity purification

For research use only.

Background:

IGFBP5, also named as IBP5, prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. It alter the interaction of IGFs with their cell surface receptors. This antibody is specific to IGFBP5.

To place an order, please [Click HERE](#).