

FTL

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 19kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human FTL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FTL; NBIA3; Ferritin light chain;

Catalog #: A1768

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2512

Isotype: IgG

Swiss Prot: P02792

Purity: Affinity purification

For research use only.

Background:

This gene encodes the light subunit of the ferritin protein. Ferritin is the major intracellular iron storage protein in prokaryotes and eukaryotes. It is composed of 24 subunits of the heavy and light ferritin chains. Variation in ferritin subunit composition may affect the rates of iron uptake and release in different tissues. A major function of ferritin is the storage of iron in a soluble and nontoxic state. Defects in this light chain ferritin gene are associated with several neurodegenerative diseases and hyperferritinemia-cataract syndrome. This gene has multiple pseudogenes.

To place an order, please [Click HERE](#).