

POLR2C

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 31kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant Protein of human POLR2C

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ei

Synonym:

RPB3; RPB31; hRPB33; hsRPB3;

Catalog #: A1785

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5432

Isotype: IgG

Swiss Prot: P19387

Purity: Affinity purification

For research use only.

Background:

This gene encodes the third largest subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. The product of this gene contains a cysteine rich region and exists as a heterodimer with another polymerase subunit, POLR2J. These two subunits form a core subassembly unit of the polymerase. A pseudogene has been identified on chromosome 21.

To place an order, please [Click HERE](#).