

GGCX

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:88kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GGCX

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VKCFD1;

Catalog #:A1806

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2677

Isotype:IgG

Swiss Prot:P38435

Purity:Affinity purification

For research use only.

Background:

This gene encodes an enzyme which catalyzes the posttranslational modification of vitamin K-dependent protein. Many of these vitamin K-dependent proteins are involved in coagulation so the function of the encoded enzyme is essential for hemostasis. Mutations in this gene are associated with vitamin K-dependent coagulation defect and PXE-like disorder with multiple coagulation factor deficiency. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).