

PTPN2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 49kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PTPN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

0

Synonym:

PTN2; PTPT; TCPTP; TC-PTP; TCELLPTP;

Catalog #: A1808

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5771

Isotype: IgG

Swiss Prot: P17706

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. Members of the PTP family share a highly conserved catalytic motif, which is essential for the catalytic activity. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. Epidermal growth factor receptor and the adaptor protein Shc were reported to be substrates of this PTP, which suggested the roles in growth factor mediated cell signaling. Multiple alternatively spliced transcript variants encoding different isoforms have been found. Two highly related but distinctly processed pseudogenes that localize to chromosomes 1 and 13, respectively, have been reported. [provided by RefSeq, May 2011]

To place an order, please [Click HERE](#).