

ICOS

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 23kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ICOS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AILIM; CD278; CVID1; MGC39850;

Catalog #: A1811

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29851

Isotype: IgG

Swiss Prot: Q9Y6W8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the CD28 and CTLA-4 cell-surface receptor family. It forms homodimers and plays an important role in cell-cell signaling, immune responses, and regulation of cell proliferation.

To place an order, please [Click HERE](#).