

SLC22A6

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 62kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SLC22A6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OAT1; PAHT; HOAT1; ROAT1;

Catalog #: A1814

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9356

Isotype: IgG

Swiss Prot: Q4U2R8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is involved in the sodium-dependent transport and excretion of organic anions, some of which are potentially toxic. The encoded protein is an integral membrane protein and may be localized to the basolateral membrane. Four transcript variants encoding four different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).