

CXADR

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 40kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CXADR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CAR; HCAR; CAR4/6;

Catalog #: A1822

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1525

Isotype: IgG

Swiss Prot: P78310

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a type I membrane receptor for group B coxsackieviruses and subgroup C adenoviruses. Several transcript variants encoding different isoforms have been found for this gene. Pseudogenes of this gene are found on chromosomes 15, 18, and 21.

To place an order, please [Click HERE](#).