

IL2RG

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 42kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human IL2RG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P64; CIDX; IMD4; CD132; SCIDX; IL-2RG; SCIDX1;

Catalog #: A1829

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3561

Isotype: IgG

Swiss Prot: P31785

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is an important signaling component of many interleukin receptors, including those of interleukin -2, -4, -7 and -21, and is thus referred to as the common gamma chain. Mutations in this gene cause X-linked severe combined immunodeficiency (XSCID), as well as X-linked combined immunodeficiency (XCID), a less severe immunodeficiency disorder. [provided by RefSeq, Mar 2010]

To place an order, please [Click HERE](#).