

BST2

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 19kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human BST2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD317; TETHERIN;

Catalog #: A1914

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 684

Isotype: IgG

Swiss Prot: Q10589

Purity: Affinity purification

For research use only.

Background:

Bone marrow stromal cells are involved in the growth and development of B-cells. The specific function of the protein encoded by the bone marrow stromal cell antigen 2 is undetermined; however, this protein may play a role in pre-B-cell growth and in rheumatoid arthritis.

To place an order, please [Click HERE](#).