

CLN3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 48kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CLN3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

BTS; JNCL;

Catalog #: A1931

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1201

Isotype: IgG

Swiss Prot: Q13286

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is involved in lysosomal function. Mutations in this, as well as other neuronal ceroid-lipofuscinosis (CLN) genes, cause neurodegenerative diseases commonly known as Batten disease or collectively known as neuronal ceroid lipofuscinoses (NCLs). Many alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).