

GLO1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 21kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GLO1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

GLYI; GLOD1

Catalog #: A1932

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2739

Isotype: IgG

Swiss Prot: Q04760

Purity: Affinity purification

For research use only.

Background:

The enzyme encoded by this gene is responsible for the catalysis and formation of S-lactoyl-glutathione from methylglyoxal condensation and reduced glutathione. Glyoxalase I is linked to HLA and is localized to 6p21.3-p21.1, between HLA and the centromere.

To place an order, please [Click HERE](#).