

RAB27A

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 24kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human RAB27A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

GS2; RAM; RAB27; HsT18676;

Catalog #: A1934

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5873

Isotype: IgG

Swiss Prot: P51159

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. The protein is membrane-bound and may be involved in protein transport and small GTPase mediated signal transduction. Mutations in this gene are associated with Griscelli syndrome type 2.

Alternative splicing occurs at this locus and four transcript variants encoding the same protein have been identified.

To place an order, please [Click HERE](#).