

EBAG9

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 24kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human EBAG9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

EB9; PDAF;

Catalog #: A1935

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9166

Isotype: IgG

Swiss Prot: O00559

Purity: Affinity purification

For research use only.

Background:

This gene was identified as an estrogen-responsive gene. Regulation of transcription by estrogen is mediated by estrogen receptor, which binds to the estrogen-responsive element found in the 5'-flanking region of this gene. The encoded protein is a tumor-associated antigen that is expressed at high frequency in a variety of cancers. Alternate splicing results in multiple transcript variants. A pseudogene of this gene has been defined on chromosome 10.

To place an order, please [Click HERE](#).