

TRIM21

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 54kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TRIM21

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1

Synonym:

RNF81; RO52; SSA; SSA1;

Catalog #: A1957

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6737

Isotype: IgG

Swiss Prot: P19474

Purity: Affinity purification

For research use only.

Background:

TRIM21 has multiple N-terminal zinc finger motifs (has E3 ligase activity), a central leucine zipper, and a potential N-glycosylation site. TRIM3 is a E3 ubiquitin-protein ligase. It can form a ubiquitin ligase complex with E2 UBE2D2, which function not only for the ubiquitination of USP4 and IKBKB but also for its self-ubiquitination. It has a role in the regulation of the cell cycle progression and could enhance the decapping activity of DCP2. It can exist in all mammalian cells as a ribonucleoprotein particle, which composed of a single polypeptide and 1 of 4 small RNA molecules.

To place an order, please [Click HERE](#).