

TFF3

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 8kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TFF3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ITF; P1B; TFI;

Catalog #: A1978

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7033

Isotype: IgG

Swiss Prot: Q07654

Purity: Affinity purification

For research use only.

Background:

Members of the trefoil family are characterized by having at least one copy of the trefoil motif, a 40-amino acid domain that contains three conserved disulfides. They are stable secretory proteins expressed in gastrointestinal mucosa. Their functions are not defined, but they may protect the mucosa from insults, stabilize the mucus layer and affect healing of the epithelium. This gene is expressed in goblet cells of the intestines and colon. This gene and two other related trefoil family member genes are found in a cluster on chromosome 21.

To place an order, please [Click HERE](#).