

SMARCA5

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP CHIP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:50 - 1:200

ChIP 1:20 - 1:100

Calculated MW: 122kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SMARCA5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

ISWI; SNF2H; hISWI; hSNF2H; WCRF135

Catalog #: A2000

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8467

Isotype: IgG

Swiss Prot: O60264

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the SWI/SNF family of proteins. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The protein encoded by this gene is a component of the chromatin remodeling and spacing factor RSF, a facilitator of the transcription of class II genes by RNA polymerase II. The encoded protein is similar in sequence to the Drosophila ISWI chromatin remodeling protein.

To place an order, please [Click HERE](#).