

BAG6

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 119kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BAG6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

G3; BAT3; BAG-6; D6S52E;

Catalog #: A2010

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7917

Isotype: IgG

Swiss Prot: P46379

Purity: Affinity purification

For research use only.

Background:

This gene was first characterized as part of a cluster of genes located within the human major histocompatibility complex class III region. This gene encodes a nuclear protein that is cleaved by caspase 3 and is implicated in the control of apoptosis. In addition, the protein forms a complex with E1A binding protein p300 and is required for the acetylation of p53 in response to DNA damage. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).