

CD83

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 23kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CD83

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BL11; HB15;

Catalog #: A2040

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9308

Isotype: IgG

Swiss Prot: Q01151

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a single-pass type I membrane protein and member of the immunoglobulin superfamily of receptors. The encoded protein may be involved in the regulation of antigen presentation. A soluble form of this protein can bind to dendritic cells and inhibit their maturation. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]

To place an order, please [Click HERE](#).