

CD69

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 20kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CD69

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AIM; EA1; MLR-3; CLEC2C; GP32/28; BL-AC/P26;

Catalog #: A2045

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 969

Isotype: IgG

Swiss Prot: Q07108

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the calcium dependent lectin superfamily of type II transmembrane receptors. Expression of the encoded protein is induced upon activation of T lymphocytes, and may play a role in proliferation. Furthermore, the protein may act to transmit signals in natural killer cells and platelets. [provided by RefSeq, Aug 2011]

To place an order, please [Click HERE](#).