

KLK3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 29kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human KLK3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

a

Synonym:

APS; KLK2A1; PSA; hK3;

Catalog #: A2052

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 354

Isotype: IgG

Swiss Prot: P07288

Purity: Affinity purification

For research use only.

Background:

Kallikrein 3 (KLK3), also known as Prostate Specific Antigen (PSA), is a member of the glandular kallikrein subfamily of serine proteases (1). It is produced by prostate epithelial cells and is secreted into prostatic ducts. Upon cleavage of 7 N-terminal amino acids (2), it is activated to liquefy semen in the seminal coagulum. Although PSA/KLK3 is produced in healthy individuals, investigators have found abnormally high levels in the blood of men with advanced prostate cancer (3,4).

To place an order, please [Click HERE](#).