

CAPZA2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 33kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CAPZA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CAPZ; CAPPA2;

Catalog #: A2054

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 830

Isotype: IgG

Swiss Prot: P47755

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the F-actin capping protein alpha subunit family. It is the alpha subunit of the barbed-end actin binding protein Cap Z. By capping the barbed end of actin filaments, Cap Z regulates the growth of the actin filaments at the barbed end. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).