

GSK3B

Reactivity: Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 48kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human GSK3B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bdimpqt

Synonym:

GSK3B; GSK3;

Catalog #: A2081

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2932

Isotype: IgG

Swiss Prot: P49841

Purity: Affinity purification

For research use only.

Background:

Glycogen synthase kinase-3 (GSK-3) was initially identified as an enzyme that regulates glycogen synthesis in response to insulin (1). GSK-3 is a ubiquitously expressed serine/threonine protein kinase that phosphorylates and inactivates glycogen synthase. GSK-3 is a critical downstream element of the PI3K/Akt cell survival pathway whose activity can be inhibited by Akt-mediated phosphorylation at Ser21 of GSK-3 and Ser9 of GSK-3 (2,3). GSK-3 has been implicated in the regulation of cell fate in Dictyostelium and is a component of the Wnt signaling pathway required for Drosophila, Xenopus, and mammalian development (4). GSK-3 has been shown to regulate cyclin D1 proteolysis and subcellular localization (5).

To place an order, please [Click HERE](#).