

MKI67

Reactivity: Human

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:200

Calculated MW: 359kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MKI67

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

KIA; MIB; MIB-1; PPP1R105;

Catalog #: A2094

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4288

Isotype: IgG

Swiss Prot: P46013

Purity: Affinity purification

For research use only.

Background:

This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X.

To place an order, please [Click HERE](#).