

CUX1

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 77kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CUX1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

μg

Synonym:

CDP; CUX; p75; CASP; CDP1; COY1; Clox; p100; p110; p200; CUTL1; GOLIM6; CDP/Cut; Cux/CDP; Nbla10317;

Catalog #: A2213

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1523

Isotype: IgG

Swiss Prot: P39880

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the homeodomain family of DNA binding proteins. It may regulate gene expression, morphogenesis, and differentiation and it may also play a role in the cell cycle progression. Several alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).