

PDCD6IP

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 96kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant Protein of human PDCD6IP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

PDCD6IP; AIP1; Alix; DRIP4; HP95; MGC17003;

Catalog #: A2215

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10015

Isotype: IgG

Swiss Prot: Q8WUM4

Purity: Affinity purification

For research use only.

Background:

Alix, a phylogenetically conserved cytosolic scaffold protein, contains an N-terminal Bro1 domain, a coiled-coil region and a C-terminal proline-rich domain (1,2). Originally identified as an ALG-2 (apoptosis-linked gene 2)-interacting protein involved in programmed cell death (3,4), Alix also regulates many other cellular processes, such as endocytic membrane trafficking and cell adhesion through interactions with ESCRT (endosomal sorting complex required for transport) proteins, endophilins, and CIN85 (Cbl-interacting protein of 85 kDa) (5,6).

To place an order, please [Click HERE](#).