

ALCAM

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:100-105kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human ALCAM

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ALCAM;CD166;MEMD;CD166 antigen;

Catalog #:A2218

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:214

Isotype:IgG

Swiss Prot:Q13740

Purity:Affinity purification

For research use only.

Background:

This gene encodes activated leukocyte cell adhesion molecule (ALCAM), also known as CD166 (cluster of differentiation 166), which is a member of a subfamily of immunoglobulin receptors with five immunoglobulin-like domains (VVC2C2C2) in the extracellular domain. This protein binds to T-cell differentiation antigene CD6, and is implicated in the processes of cell adhesion and migration. Multiple alternatively spliced transcript variants encoding different isoforms have been found. [provided by RefSeq, Aug 2011]

To place an order, please [Click HERE](#).