

CHD3

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 227kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human CHD3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A2221

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1107

Isotype: IgG

Swiss Prot: Q12873

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).