

TRIM28

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:3000

Calculated MW: 89kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TRIM28

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fl

Synonym:

KAP1; TF1B; RNF96; TIF1B;

Catalog #: A2245

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10155

Isotype: IgG

Swiss Prot: Q13263

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene mediates transcriptional control by interaction with the Kruppel-associated box repression domain found in many transcription factors. The protein localizes to the nucleus and is thought to associate with specific chromatin regions. The protein is a member of the tripartite motif family. This tripartite motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region.

To place an order, please [Click HERE](#).