

## IDH1

---

**Reactivity:** Human Mouse

**Tested applications:** WB IHC IF FC

**Calculated MW:** 414

**Observed MW:** Refer to figures

**Immunogen:**

A synthetic Peptide of human IDH1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:** A2263

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 3417

**Isotype:** IgG

**Swiss Prot:** O75874

**Purity:** Affinity purification

For research use only.